

Free Stand – Product Cleaning and Care

- The Free Stand work surface is made of ABS plastic blend and the work surface is coated with soft-touch paint. The frame is made of extruded aluminum and powder-coat painted.
- Exercise caution with the work surface in proximity to sharp objects and when setting a unit down on an abrasive surface such as pavement or concrete. Contact with sharp objects or abrasive surfaces can cause scuffing and scratching of the work surface.
- Clean the aluminum frame and work surface frequently by applying a household cleaner (such as Formula 409® or Green Works®) and wiping dry with a clean, soft cloth. Do not allow cleaner chemicals to remain on the work surface for an extended period of time.

Recommended Cleaners

- Formula 409
- Green Works
- Windex
- Bleach and Water solution (10% diluted)
- HB Quat
- TB Quat
- Virex II
- Lysol Disinfecting Wipes

Avoid Using – can cause mild staining and discoloration

- Clorox Disinfecting Wipes
- Sanizide

Never Use – can cause severe staining, discoloration, or peeling of paint

- Hand Sanitizer
- Mr. Clean Magic Eraser
- Ethyl Alcohol